


Recovery time after opening the door

At the highest test temperature and when the temperature has stabilised, the door of the oven is opened for 30 seconds and then the time is measured until the **actual temperature** in the middle of the working space remains permanently within the limit deviation specified in DIN 12880:2007-05.



End of the recovery time (DIN 12880:2007-05)

[Overview Glossary Temperature control chamber](#)

Autor:

www.atmosafe.net > [Glossary](#) > [DIN 12880:2007-05](#) > [Recovery time](#)

AtmoSAFE is a brand of Memmert GmbH + Co. KG
Copyright © 2009 Memmert GmbH + Co. KG.
All Rights Reserved.



memmert
Experts in Thermostatics