

# Setting and programming parameters

**Digital time switches have replaced mechanical timers in the modern temperature control chamber.**

Apart from the manual input of a specific programme runtime, repetitions, specific weekdays, temperature ramps (time-dependent temperature progressions) as well as different combinations of these options, can be entered directly on the appliance, or via software, depending on the type of appliance. Apart from temperature, further parameters such as fan speed, air exchange rate, pressure, humidity or CO2 content can be digitally controlled in many appliances, and thus integrated within the ramp programming.

[Overview](#) [Glossary](#) [Temperature control chamber](#)

Picture credit: Memmert GmbH + Co. KG

Autor:

---

[www.atmosafe.net](http://www.atmosafe.net) > [Glossary](#) > [Operation](#) > [Programming](#)

---

AtmoSAFE is a brand of Memmert GmbH + Co. KG  
Copyright © 2009 Memmert GmbH + Co. KG.  
All Rights Reserved.



Memmert humidity chamber HCP



**memmert**  
Experts in Thermostatics